

## Product Description

### Vista Patinas and Traditional Patinas Sculpt Nouveau

Vista Patinas are traditional patinas because they contain chemicals which react with the metal. They contain a base blend of chemicals that produce a beautiful green on bronze, brass, and copper, or on iron & steel they will create exquisite rust.

Added to this traditional base solution are different colored oxides, which will blend with the chemical reaction to form a marbled patina. For example, using the Vista Black on bronze would produce a green and black patina – on iron there would be rust and black together. These patinas are not recommended for aluminum. The color will vary with different metals.

#### Vista Patinas on Bronze, Brass, or Copper



**Vista Blue**



**Vista Green**



**Vista Yellow**



**Vista Rust**



**Vista Burgundy**



**Vista Red**



**Vista Brown**



**Vista Black**

#### Vista Patinas on Iron & Steel

